

# Endoscopic reprocessing: TOWARDS STERILIZATION?

## SPAULDING CLASSIFICATION IN MODERN ENDOSCOPY

As per CDC definition<sup>1</sup> Spaulding Classification is “a strategy for sterilization or disinfection of inanimate objects and surfaces based on the degree of risk involved in their use.” There are three groupings of medical devices: lower risk non-critical medical devices which make contact with outer body surface, e.g. intact skin and require low to intermediate level disinfection for safe use; semi-critical devices, which make contact with broken skin and inner body surface, e.g. mucous membranes lining body cavities, these devices must be at least high level disinfected and preferably sterilized after each use; critical devices which make contact with blood, organs and sterile body tissues, these devices must be sterile for each use.



SPEAKER:

## Olegs Rasheed Tucs

EDUCATOR, CSSD MANAGER AND RESEARCHER WITH WIDE EXPERIENCE IN HEALTHCARE ENVIRONMENTS

MASTER OF SCIENCES IN CELL AND MOLECULAR BIOLOGY

BACHELOR OF SCIENCES IN DECONTAMINATION MANAGEMENT

CERTIFIED SENIOR EMBRYOLOGIST (ESHRE)



Call Out#  
“Different procedures  
have varied levels of  
infection transmission  
risk.”

RISK ASSESSMENT  
SHOULD BE DONE FOR  
ENDOSCOPIC PROCEDURES  
TO CHOOSE APPROPRIATE  
REPROCESSING STRATEGY.



## CRITICAL ENDOSCOPIC PROCEDURES

Rigid endoscopes are used in laparoscopic and arthroscopic surgical procedures and are classified in the same critical risk category as invasive surgical instruments. Semi-flexible and flexible endoscopes are routinely used for interventional invasive medical procedures such as biopsies, polypectomies, injection etc. accessing sterile body tissues and circulatory systems by breaching mucous membranes. Single use sterile tools such as snares, brushes, needles, and others are passed through endoscope biopsy channels before coming in contact with sterile body tissues.



**Call Out#**  
**“Invasive procedures can be performed by accessing sterile body tissues through mucous membranes.”**

ENDOSCOPES ARE ROUTINELY USED FOR CRITICAL INVASIVE PROCEDURES.

## DIFFERENCES BETWEEN THE HLD AND STERILIZATION

Scopes can carry as much as  $10^7$ - $10^{10}$  bacteria after the procedure<sup>2</sup>. Depending on risk classification endoscopes must be high level disinfected or sterilized for each use. High-level disinfection complete elimination of all microorganisms in or on an instrument, except for small numbers of bacterial spores. The CDC defines high-level disinfectant as a sterilant used for a shorter contact time to achieve a  $6\text{-log}_{10}$  kill of an appropriate Mycobacterium species<sup>3</sup>. Sterilization – killing of all viable pathogens. The chance of a single viable pathogen surviving sterilisation process is less than one in a million, e.g. the Sterility Assurance Level is  $10^{-6}$



## Call Out#

**“High-level disinfection significantly reduces endoscope contamination but there can be surviving microorganisms post disinfection including bacterial spores.”**



STERILIZATION IS  
VALIDATED TO  
KILL ALL PATHOGENS  
ON AN INSTRUMENT.

## THE BEST PRACTICE FOR PATIENT SAFETY

There are numerous studies proving survival of pathogenic microorganisms on endoscopes after high level disinfection. These disinfection survivors have caused numerous outbreaks of healthcare associated infections and led to procedure complications including fatalities<sup>4,5,6</sup>. Meticulous cleaning followed by sterilization or use of sterile single use endoscopes can substantially reduce risk of infection transmission and improve patient safety.



## Call Out#

**“Including sterilization as part of endoscope reprocessing where possible can eliminate residual contamination and make endoscopic procedures safer for patients.”**

BEST PRACTICE MAY  
MEAN TO GO ABOVE THE  
LOWEST ACCEPTABLE  
DECONTAMINATION  
LEVEL.

# Take home messages:

- ▶ ENDOSCOPES ARE ROUTINELY USED FOR INTERVENTIONAL INVASIVE MEDICAL PROCEDURES BRINGING THEM IN CRITICAL HIGH-RISK CATEGORY.
- ▶ BACTERIA CAN SURVIVE IN ENDOSCOPES AND CAUSE INFECTION OUTBREAKS DESPITE CAREFUL REPROCESSING OF INSTRUMENTS IN ACCORDANCE WITH MANUFACTURER'S GUIDELINES.
- ▶ HIGH-LEVEL DISINFECTION SIGNIFICANTLY REDUCES ENDOSCOPE CONTAMINATION BUT THERE CAN BE SURVIVING MICROORGANISMS POST DISINFECTION INCLUDING BACTERIAL SPORES.
- ▶ EVEN THOUGH HIGH-LEVEL DISINFECTION MAY BE ACCEPTABLE METHOD OF ENDOSCOPE REPROCESSING STERILIZATION WHERE POSSIBLE CAN ELIMINATE RESIDUAL CONTAMINATION AND MAKE ENDOSCOPIC PROCEDURES SAFER FOR PATIENTS.

Bronchoscope image: <https://www.olympus-europa.com/medical/en/Products-and-Solutions/Products/Product/BF-XT190.html>

References: **1.** <https://www.cdc.gov/infectioncontrol/guidelines/disinfection/rational-approach.html> **2.** Rutala WA, Weber DJ, 2016. Outbreaks of carbapenem resistant Enterobacteriaceae infections associated with duodenoscopes: what can we do to prevent infections?. Am J Infect Control. 2016 May 2;44(5 Suppl):e47-51. doi: 10.1016/j.ajic.2015.10.037. **3.** <https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines-H.pdf> **4.** Epstein L, Hunter JC, Arwady MA, Tsai V, Stein L, Gribogiannis M, Frias M, Guh AY, Laufer AS, Black S, Pacilli M, Moulton-Meissner H, Rasheed JK, Avellan JJ, Kitchel B, Limbago BM, MacCannell D, Lonsway D, Noble-Wang J, Conway J, Conover C, Vernon M, Kallen AJ. New Delhi metallo-β-lactamase-producing carbapenem-resistant Escherichia coli associated with exposure to duodenoscopes. JAMA. 2014 Oct 8;312(14):1447-55. doi: 10.1001/jama.2014.12720. PMID: 25291580. **5.** Ofstead CL, Quick MR, Wetzler HP, Eiland JE, Heymann OL, Sonetti DA, Ferguson JS. Effectiveness of Reprocessing for Flexible Bronchoscopes and Endobronchial Ultrasound Bronchoscopes. Chest. 2018 Nov;154(5):1024-1034. doi: 10.1016/j.chest.2018.04.045. Epub 2018 May 31. PMID: 29859183. **6.** <https://www.fda.gov/medical-devices/letters-health-care-providers/infections-associated-reprocessed-urological-endoscopes-letter-health-care-providers>.

**ASP** Advanced Sterilization Products

ASP International GmbH, Zug Branch  
Bahnhofstrasse 2, Zug 6300, Switzerland  
©ASP 2022. All Rights Reserved.

[asp.com](http://asp.com)



ADVANCED STERILIZATION PRODUCTS, INC.  
33 Technology Drive, Irvine CA 92618, USA



ASP, The Netherlands BV  
BIC 1, 5657 BX Eindhoven, The Netherlands



0123